

## ADHERENCE ASSAY (revised 8/23/02)

### Preparing Cells

\*\*prepare a 24 well plate with sterile glass coverslips

- 1) Grow A549 cells to ~80-90% confluency in a T-75
- 2) Trypsinize cells, spin down and remove media by aspiration
- 3) Resuspend cells in 5 ml complete growth media
- 4) Mix 10  $\mu$ l of cell suspension with with an equal volume of trypan blue
- 5) Add 10  $\mu$ l to hemocytometer and count cells in four corners
- 6) Average four counts to calculate the cell concentration:

$$2 \times \text{average cell count} \times 10^4 = \text{number of cells/ml}$$

- 7) Calculate the volume of cells required to give  $2 \times 10^5$  cells/well

volume of cells required (ml) =  $(2 \times 10^5 \text{ cells})(\# \text{wells})/\text{conc from hemocytometer}$

- 8) Adjust volume to give  $2 \times 10^5$  cells/ml (ie 1 ml per well)
- 9) Add 1 ml of A549 cells into each well with coverslip
- 10) Allow cells to adhere overnight

### Binding Assay

- 1) Remove media and wash A549 cells 2X with HBSS
- 2) Add 1 ml of complete growth media WITHOUT antibiotics or FBS to each well and incubate at 37 °C, 5% CO<sub>2</sub>
- 3) Grow bacterial culture to OD<sub>600</sub> ~ 1.0

Pseudomonas: An OD<sub>600</sub> of 1.0 is equal to  $5 \times 10^8$  cfu/ml  
spin down cells and resuspend in cell culture media to give a final concentration of  $10^7$  cfu/ml

- 4) Remove media from each well
- 5) Add 1 ml of bacterial suspension to each well
- 6) Incubate at 37 °C for desired amount of time
- 6) Fix coverslip and stain